CELLBOOSTER® LIFT:
Hyaluronic Acid 18 mg Rejuvenating Complex

INSTRUCTIONS FOR USE

DESCRIPTION

CELLBOOSTER® LIFT: Hyaluronic Acid Rejuvenating Complex is produced using patented CHAC Technology. It consists of non-crosslinked Hyaluronic Acid with Amino acids and Vitamins. The device is designed to diminish the effects of aging by filling the soft tissue, preventing oxidative stress and providing a deep hydration. It stimulates collagen production, improves microcirculation and cell nutrition.

The device is a clear, yellow, sterile, injectable solution for single use presented in a case of 6 x 3 mL vials.

COMPOSITION

Each vial contains:
 • non-crosslinked Sodium Hyaluronate 18 mg;
 • amino acids (Arginine, Glycine, Lysine, Proline, Valine);
 • vitamins (Riboflavin, Biotin, Sodium Ascorbyl Phosphate, Tocopherol);
 • phosphate buffer pH 6.8-7.6 q.s.

INDICATIONS

The device is a resorbable injectable implant for use in the epidermis or dermis for diminishing moderate skin depression and providing dermis re-densification and provides deep hydration, and improvement with microcirculation. Reduces dryness and strengthens the extracellular matrix. Further it diminishes the effects of aging by filling the soft tissue and preventing oxidative stress. Restores the structure of the skin and subcutaneous tissue. Evens out skin texture. Effective for the attenuation of wrinkles and skin laxity, for a lifting effect and skin tone improvement. Restores skin with atrophic scars, post-acne, striae. Prepares skin for plastic surgery

CONTRAINDICATIONS

The device must not be used in patients:
 • with acute infectious diseases;
 • with severe somatic and mental illnesses;
 • with autoimmune and allergic diseases;
 • with known allergic or hypersensitive reactions to the components of the product;
 • with varices in the area;
 • with inflammations or infections in the area of the treatment;
 • who tend to develop keloids and hypertrophic scarring;
 • who received temporary or permanent implants in the same area;
 • who receives anticoagulant treatment;
 • who simultaneously receives laser treatment, deep chemical peels or dermabrasion;
 • pregnant or breastfeeding women;
 • children.

PRECAUTIONS FOR USE

The device must only be used by physicians who have received special training in the mesotherapy injection techniques. The device must be injected into non-inflamed, disinfected, healthy skin. All asepsis and hygiene precautions must be taken to prepare the site, equipment and patient before implantation. Sensitive skin may be pretreated using a local anaesthetic patch or cream. After use, discard the syringe, needle and remaining product in a special container.

UNDESIRABLE SIDE EFFECTS

The physician must inform the patient about potential side effects associated with implantation of this device that may occur immediately or within a reasonable delay.

The following events and reactions have been observed with the device or similar products:
 • stinging, redness or slight local inflammation which should disappear after 24 to 48 hours;
 • mild edema and small hematomas may develop but should disappear within 48 hours.

Patients must inform their physician as soon as possible about any other undesirable effect, or the effects listed above that last longer than the specified duration. The physician should treat these side effects appropriately. Such undesirable side effects must be reported to the distributor and/or manufacturer.

METHOD OF USE

The device must be injected in the epidermis or dermis of the face, neck, decollete area, back of the hands, internal face of the arms or body according to the mesotherapy technique.

The device must be implanted aseptically, using single-use sterile syringe which is not supplied by Suisselle.

Duration of treatment:
 • intensive course: 3 treatments with 2-3 weeks intervals;
 • support course: 1 treatment per month.

WARNINGS

• Before use, verify the expiry date and integrity of the package; do not use after the expiry date or if the package is damaged.
 • Do not use if the solution is cloudy or its color differs from that described above.
 • Do not mix with other products.
 • Do not inject into blood vessels, bones, tendons, ligaments, nerves or muscles.

- Do not overcorrect.
- Once the vial has been opened, use immediately.
- Vial is intended for a single use: do not reuse and do not re-sterilize. Quality and sterility can only be guaranteed for a vial in its original packaging. The re-use of the device creates a potential infection risk for patients or users.
- Patients should be advised not to apply any make-up for 12 hours after the injection and to avoid prolonged exposure to the sun, UV rays, frost, as well as using sauna or hammam, for two weeks following implantation.
- There are incompatibilities between Sodium Hyaluronate and Quaternary Ammonium Compounds such as Benzalkonium Chloride. Therefore, the device should never be placed into contact with these substances, or with medical-surgical instruments that have been treated with these substances.

STORAGE CONDITIONS

The device should be stored at 5-25°C in a dry place in the original box and protected from direct sunlight, heat or frost. Handle with care.

EXPLANATION OF INTERNATIONAL SYMBOLS

| | |
|--|---|
| | Consult instructions for use |
| | Use by |
| | Do not use if package is damaged |
| | Temperature limitation |
| | Do not re-sterilize |
| | Keep dry |
| | Handle with care |
| | Medical Device |
| | Single sterile barrier system with protective packaging outside |
| | Authorized representative in the European Community |
| | Email: info@suisselle.com |
| | Made in Switzerland |

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